# Column for reversed phase

# C30 Series

# Develosil C30-UG [Develosil RPAQUEOUS, RPFULLERENE, Combi-RP] Develosil RPAQUEOUS-AR

It is chemical modification about a triacontyl radical (c30) based on a high grade silica gel.

# Physical properties of Develosil C30Series

Column name	Ligand	Carbon	End Capping	Surface area	Pore diameter	Pore volume	Range of pH
C30-UG							
RPAQUEOUS	Triacontyl radical (Monfunctional)		Yes (Double) 300m	300m²/g	lm²/g 14nm	1.05mL/ g	pH2-8
RPFULLERENE		18%					
Combi-RP		10 /0		300iii /g			
RPAQUEOUS-AR	Triacontyl radical (Trifunctional)						pH1-7


**※**C30-UG, RPAQUEOUS, RPFULLERENE, and Combi-RP are using the same Packing material. The name of article is changed and shipped to purpose-oriented.

#### The feature of Develosil C30 series

Different characteristics from an ODS column	To solution well inseparable using ODS.  Develosil C30 series does not only simply have a long alkyl chain.
A 100% of water mobile phase can also be used	It is the present condition for retention not to become early rapidly in the usual column, if a 100% of water mobile phase is used, or not to be stabilized.  In Develosil C30 series, by controlling a base material surface well shows high reproducibility also by a 100% of water mobile phase.
Solid recognition ability is high	Solid recognition ability changes with a monomeric type or a polymeric type, and the length of an alkyl chain is also one of factors.  It is fit for an isomer, separation of polycyclic aromatic, etc.
The best for separation of a liposolubility compound	A result good also for carotenoid (especially carotene) or fullerene can be obtained.

## The difference between Develosil C30-UG (RPAQUEOUS) and RPAQUEOUS-AR

#### Comparison 1 of column performance




Column	: Develosil® C30-UG-5
	Develosil® RPAQUEOUS-AR-5
Size	: 4.6x150mm
Mobile phase	: MeOH/Water=70/30
Flow rate	: 1.0ml/min
Temperature	: 40°C
Detection	: UV254nm
Sample	: 1.Uracil / 2.Caffeine / 3.Phenol
	4.Methyl Benzoate / 5.Benzene
	6.Toluene

The left figure is data based on the packingmaterial appraisal method of our regulation. From this chromatogram, a hydrogen associativity and hydrophobic surface polarity are evaluated.

This result shows that the equivalent degree of separation is obtained also by C30-UG or RPAQUEOUS-AR.

ICOIUMN name	Hydrogen associativity	Hydrophobicity	Surface polarity
C30-UG	0.47	1.73	0.44
RPAQUEOUS-AR	0.48	1.74	0.42

#### Comparison 2 of column performance



Conditions;	
Column	: Develosil® C30-UG-5
	Develosil® RPAQUEOUS-AR-5
Size	: 4.6x150mm
Mobile phase	: MeOH/Water=80/20
Flow rate	: 1.0ml/min
Temperature	: 40°C
Detection	: UV254nm
Sample	: 1.Uracil / 2.o-Terphenyl / 3.Triphenylene

It is comparison of solid recognition ability. It will be said that the power of recognizing a solid compound is strong, so that this figure is high.


C30-UG-5 α=1.58

RPAQUEOUS-AR α=1.73

It turns out that the direction of RPAQUEOUS-AR is in the tendency to recognize a solid compound.

#### Durability of the column in 100% of water

It is a data figure in an acid-proof examination. RPAQUEOUS-AR which has adopted polymeric type C30 has high durability on acid terms. An ion pair reagent can be used and a result good for an amino acid, analysis of peptide, etc. can be obtained.




Mobile phase: 0.5%TFA Temperature: 40°C

#### Durability of the column at the case of a contain organic solvent

The data of durability test at the time of stepping on an organic solvent to a mobile phase-proof is shown.

Under these terms, use of nourishment is enabled by the acid-proof examination for 500 hours.




Mobile phase: CAN / 0.1 % Phosphoric acid = 98/2

Temperature: 40°C

#### Durability under 100% of a basin system, and alkali conditions

The data in an alkali resistance test is shown. High durability is shown under the alkali terms of 100% of a basin system.


Deterioration condition differs in the mobile phase to be used with pH and temperature.



# **Application**


# Develosil RPAQUEOUS, RPAQUEOUS-AR

# Amino acid I



Column	Develosil RPAQUEC		US
		4.6x250mr	n
Mobile phase	Water		
Flow rate	1.0mL/min		
Temperature	40°C		
Detection	RI		
Infection rate	10µL		
Sanmple	1=Alanine	0.2%	
	2=Valine	0.2%	
	3=Isoleuci	ne 0.2%	
	4=Leucine	0.2%	


#### Amino acid II



	Column	Develosil RPAQUEOUS-AR-5
		4.6x150mm
	Mobile phase	0.1% TFA or 0.1% HFBA
	Flow rate	1.0mL/min
	Temperature	30°C
	Detection	UV@210nm or RI
	Sanmple	1=D-Alanine
		2=D-Alanine-D-Alanine


The capability of an HFBA of an ion pear is high and it also comes to hold alanine.

### Catecholamine




Column	Develosil RPAQUEOUS
	4.6x250mm
Mobile phase	0.1% M KH2PO4,
	pH2.6 with H3PO4
Flow rate	1.0mL/min
Temperature	30°C
Detection	UV@210nm
Sanmple	1=Norepinephrine
	2=Epinephrine
	3=Dopamine
	4=DL - Dopa

### **Nucleotide**




Column Develosil RPAQUEOUS 4.6x250mm Mobile phase 0.1% M KH2PO4, pH6.0 with KOH 1.0mL/min Flow rate 30°C Temperature UV@260nm Detection Sanmple 1=5'-CTP 2=5'-CDP 3=5'-CMP 4=5'-GTP 5=5'-GDP 6=5'-GMP 7=5'-ATP 8=5'-ADP 9=5'-AMP

# Water soluble vitamin




Column	Develosil RPAQUEOUS
	4.6x150mm
Mobile phase	50M Na <sub>3</sub> PO <sub>4</sub> (pH6.9)
Flow rate	1.0mL/min
Temperature	30°C
Detection	UV@254nm
Sanmple	1=Orotic acid
	2=Nicotinic acid
	3=Pyridoxal
	4=Pyridoxine
	5=Nicotinamide

# Organic acid




Column	Develosil RPAQUEOUS-AR-5	
	4.6x150mm	
Mobile phase	0.1% TFA	
Flow rate	1.0mL/min	
Temperature	30°C	
Detection	UV@210nm	
Sanmple	1=Pyruvic acid (0.003%)	
	2=Lactic acid (0.06%)	
	3=Citric acid (0.08%)	

# The organic acid in a kiwi berry



	Column	Develosil RPAQUEOUS-AR-3
		4.6x250mm
	Mobile phase	0.1% Phosphoric acid
	Flow rate	0.6mL/min
	Temperature	20°C
)	Detection	UV@210nm
	Sanmple	1=Quinic acid
		2=Malic acid
		3=Citric acid

## Ascorbic acid




Column	Develosil RPAQUEOUS-AR-5
	4.6x250mm
Mobile phase	0.1% TFA
	(Trifluoroacetic acid)
Flow rate	1.0mL/min
Temperature	30°C
Detection	UV@254nm
Sanmple	1=Ascorbic acid


C30 stationary phase holds ascorbic acid.
Retention (k) is set to about 0.8 because t0 is 1.9 minutes.

# **Application**

# Develosil RPAQUEOUS, RPAQUEOUS-AR




## Develosil C30-UG



#### Develosil Combi-RP

### Separation of drinkable preparations


(High throughput analysis)



Column	Develosil Combi-RP 4.6x50mm
	Develosil ODS 4.6x50mm
Mobile phase	A) 0.1% TFA
	B) Acetonitrile
	% B 0-100% 5 minutes
Flow rate	3.0mL/min
Temperature	30°C
Detection	UV@254nm

#### Develosil RPFULLERENE

#### The diastereomer of sugar addition fullerene



Column	Develosil RPFULLERENE
	4.6x250mm
Mobile phase	Toluene/acetonitrile = 4:6
Flow rate	1.0mL/min
Detection	UV@335nm