The packed column for HPLC

Develosil XG-CN

The 3rd phase "XG-CN" of Develosil XG series Improvement point: The durability over acetonitrile was improved sharply.

The trait of Develoseil XG-CN

Silica base

This product modify the cyanopropyl group based on a high grade silica gel.

The column which it is easy to use As for this column, equilibrium is early and the pressure is low, too.

In addition, even a normal phase is available in the reverse phase.

Characteristic of Develosil XG-CN

Surface area	300 m ² /g
Pore diameter	14nm
Pore volume	1.10 mL/g
Ligand	Cyanopropyl
Carbon content	7.50%
Endcapping	0
pH range	pH 1.5 - 8

Stability

Develosil keeps high stability to undertake all processes in own company.


Use solvent

Compared with the conventional CN column, the durability over acetonitrile improves considerably.

In addition, I do not choose the kind of the solvent.

2. Performance of Develosil XG-CN

Organic solvent composition and a pressure




Develosil XG series shows low pressure conventionally. This became an ease of use column. The flow rate of the above terms can be raised and used. And, a solvent with high viscosity can also be used.

Use by acetonitrile is possible

In the conventional CN column, the specification of acetonitrile had restriction. In Develosil XG-CN, the durability over acetonitrile is improving greatly. I can thereby use it without a problem.


Difference



The MeOH system was used for the performance test at the time of shipment until now. But in Develosil XG-CN, the performance test is done with the solvent of an ACN system.

3. Develosil XG-CN performance evaluation

Performance evaluation of time of shipment.




Conditions;
Column : Develosif XG-CN (5um)
Size : 4.6x150mm
Mobile Phase : ACN/Water=50/50

Flow rate : 1.0ml/min
Temperature : 30°C
Detection : UV254nm
Sample : 1. Benzene
2. 1,3,5 Triphenyl benzene

Time of shipment performance test is evaluated in an ACN system. Only the column which passed the basis is shipped.


Hydrogen bonding capacity • Hydorphobicity • Surface polarity



Conditions;
Column : Develosit® XG-CN (5um)
Size : 150x4.6mm
Mobile Phase : CH₃OH/Water (30:70)
Flow rate : 1.0ml/min
Temperature : 40°C
Detection : UV 254nm
Sample : 1. Uracit/2. Caffeine
3. Phenol/4. Benzene
5. Methyl benzoate/6. Toluene

Hydrogen bonding capacity α "Caffeine/Phenol", hydrophobicity α "Toluene/Benzene", and surface polarity α "Methyl Benzoate/Toluene) are estimation.

Plane selectivity



Conditions;

Column : Develosil® XG-CN (5um)
Size : 4.6x150mm

Mobile Phase: MeOH/Water=65/35
Flow rate: 1.0ml/min
Temperature: 40°C


Detection : UV254nm

Sample : 1. Uracil / 2. o-Terphenyl

Triphenylene

Plane selectivity α "Triphenylene/o-Terphenyl" is estimation.

Coordinate compounds




Conditions;
Column : Develosil XG-C18M (5um)
Size : 150x4.6mm
Mobile Phase : CH₃CN/0.2% H₃PO₄ (35:65)
Flow rate : 1.0ml/min
Temperature : 40°C
Detection : UV254nm

Sample : 1. Uracil / 2. Toluene
3. Quinizarine

Plane selectivity α "Quinizarine/Toluene" is estimation.

4. Application data

Analysis of steroid



Conditions;

Column : Develosit® XG-CN (5um)

Size : 4.6x150mm

Mobile Phase : ACN/Water=35/65

Flow rate : 1.0ml/min


Temperature : 40°C

Detection : UV210nm Sample : 1. Estriol (0.13mg/ml)

β-Estradiol (0.13mg/ml)
 17α-Estradiol (0.13mg/ml)

Estrone (0.13mg/ml)
 Injection volume : 3.0uL

Analysis of a base compound



Column : Develosil XG-C18M (5um)

ize : 150x4.6mm


Mobile Phase : CH₃CN/50mM Ammonium Acetate (pH7.0) (20/80)

Flow rate : 1.0ml/min
Temperature : 40°C

Detection : UV254nm Sample : 1. Uracii (0.08mg/ml)

2. Pyridine (0.13mg/ml) 3. Phenol (0.21mg/ml)

Analysis of Acidic compunds



Conditions;

Column : Develosil® XG-CN (5um)
Size : 4.6x150mm
Mobile Phase : ACN/0.2%

Phosphoric acid=2/98


Flow rate : 1.0ml/min Temperature : 40°C

Detection : UV254nm
Sample : 1. Ascorbic acid (0.02mg/ml)
2. Benzoic acid (0.16mg/ml)

Acetylsalicylic acid (0.16mg/ml)

Salicylic acid (0.16mg/ml)
Injection volume : 2.0uL

The example of analysis using normal phases system solvent



Conditions;
Column : Develosil® XG-CN (5um)

Size : 4.6x150mm

Mobile Phase : Heptane
Flow rate : 1.0ml/min
Temperature : 30°C

Detection : UV254nm

Sample : 1. Anthracene (0.03mg/ml)

2. Pyrene (0.51mg/ml) 3. Chrysene (0.26mg/ml)

4. Benzo [a] Pyrene (0.26mg/ml) Injection volume: 10uL

Use by a normal phases can be performed also in Develosil XG-CN. Specification of a "reversed phase" and a "normal phases" is possible at the time of an order.

5. Price List

I.D.2.0mm,3.0mm,4.0mm same price

Convention columns

■ Conventiona columns 【Reversed phase】

Packing name	Column Size	Particle size	P/N	Price (yen)
	I.D 4.6x35mm		XGCN346035W	38,000
	I.D 4.6x50mm		XGCN346050W	38,000
	I.D 4.6x75mm	3	XGCN346075W	43,000
	I.D 4.6x100mm	၁	XGCN346100W	48,000
	I.D 4.6x150mm		XGCN346150W	58,000
XG-CN	I.D 4.6x250mm		XGCN346250W	68,000
AG-CIV	I.D 4.6x35mm	5	XGCN546035W	35,000
	I.D 4.6x50mm		XGCN546050W	35,000
	I.D 4.6x75mm		XGCN546075W	40,000
	I.D 4.6x100mm		XGCN546100W	45,000
	I.D 4.6x150mm		XGCN546150W	55,000
	I.D 4.6x250mm		XGCN546250W	65,000

Guard column for conventiona

[Reversed phase]

Packing name	P/N	Price (yen)
Gurad cartridge(4piceces) 4.0x10mm	XGCN540010C	18,000
Holder 4.0x10mm	HO00040010C	9,500
Guard cartridge set ※	XGCN540010W	14,000
☆ Guard cartridge set = Holder (1piece) + cartridge 4.0x10m	nm (1piece)	

◆ Conventiona columns 【Normal phase】

Packing name	Column Size	Particle size	P/N	Price (yen)
	I.D. 4.6x35mm		XGCN346035WN	38,000
	I.D. 4.6x50mm		XGCN346050WN	38,000
	I.D. 4.6x75mm	3	XGCN346075WN	43,000
	I.D. 4.6x100mm	ა	XGCN346100WN	48,000
	I.D. 4.6x150mm		XGCN346150WN	58,000
XG-CN	I.D. 4.6x250mm		XGCN346250WN	68,000
AG-CIV	I.D. 4.6x35mm		XGCN546035WN	35,000
	I.D. 4.6x50mm		XGCN546050WN	35,000
	I.D. 4.6x75mm	_	XGCN546075WN	40,000
	I.D. 4.6x100mm	5	XGCN546100WN	45,000
	I.D. 4.6x150mm		XGCN546150WN	55,000
	I.D. 4.6x250mm		XGCN546250WN	65,000

Guard column for conventiona

[Normal phase]

Packing name	P/N	Price (yen)
Gurad cartridge(4piceces) 4.0x10mm	XGCN540010CN	18,000
Holder 4.0x10mm	HO00040010C	9,500
Guard cartridge set ※	XGCN540010WN	14,000
	ım (1piece)	

Semi-micro columns

■ Semi-micro columns 【Reversed phase】

Packing name	Column Size	Particle size	P/N	Price (yen)
	I.D. 2.0x35mm		XGCN320035W	38,000
	I.D. 2.0x50mm		XGCN320050W	38,000
	I.D. 2.0x75mm	3	XGCN320075W	43,000
	I.D. 2.0x100mm		XGCN320100W	48,000
	I.D. 2.0x150mm		XGCN320150W	58,000
XG-CN	I.D. 2.0x250mm		XGCN320250W	68,000
VG-CIA	I.D. 2.0x35mm		XGCN520035W	35,000
	I.D. 2.0x50mm		XGCN520050W	35,000
	I.D. 2.0x75mm	_	XGCN520075W	40,000
	I.D. 2.0x100mm	5	XGCN520100W	45,000
	I.D. 2.0x150mm		XGCN520150W	55,000
	I.D. 2.0x250mm		XGCN520250W	65,000

Guard column for Semi-micro

[Reversed phase]

Packing name	P/N	Price (yen)
Gurad cartridge(4piceces) 1.5x10mm	XGCN515010CN	20,000
Holder 1.5x10mm	HO00015010C	15,000
Guard cartridge set ※	XGCN515010WN	20,000
※ Guard cartridge set = Holder (1piece) + cartridge 1.5x 10m	nm (1piece)	

◆ Semi-micro columns 【Normal phase】

Packing name	Column Size	Particle size	P/N	Price (yen)
	I.D. 2.0x35mm		XGCN320035WN	38,000
	I.D. 2.0x50mm		XGCN320050WN	38,000
	I.D. 2.0x75mm	3	XGCN320075WN	43,000
	I.D. 2.0x100mm	3	XGCN320100WN	48,000
	I.D. 2.0x150mm		XGCN320150WN	58,000
XG-CN	I.D. 2.0x250mm		XGCN320250WN	68,000
VG-CIA	I.D. 2.0x35mm		XGCN520035WN	35,000
	I.D. 2.0x50mm		XGCN520050WN	35,000
	I.D. 2.0x75mm	_	XGCN520075WN	40,000
	I.D. 2.0x100mm	5	XGCN520100WN	45,000
	I.D. 2.0x150mm		XGCN520150WN	55,000
	I.D. 2.0x250mm		XGCN520250WN	65,000

Guard column for Semi-micro

[Normal phase]

Packing name	P/N	Price (yen)
Gurad cartridge(4piceces) 1.5x10mm	XGCN515010C	20,000
Holder 1.5x10mm	HO00015010C	15,000
Guard cartridge set 🔆	XGCN515010W	20,000
	nm (1piece)	

5. Price List

Prep columns

■ Prep columns 【Reversed phase】

Packing name	Column Size	Particle size	P/N	Price (yen)
XG-CN	I.D. 10.0x150mm	5	XGCN5P1150W	130,000
	I.D. 10.0x250mm		XGCN5P1250W	150,000
	I.D. 20.0x150mm		XGCN5P2150W	240,000
	I.D. 20.0x250mm		XGCN5P2250W	275,000

Packing name	Column Size	P/N	Price (yen)
	I.D. 8.0x10mm	XGCN580010W	20,000
XG-CN	(for I.D. 10.0mm)	AGCINSOUD TOVV	
	I.D. 20.0x50mm XGCN5P2050W		105,000
The above goods			
It becomes a packing guard column.			

◆ Prep columns 【Normal phase】

Packing name	Column Size	Particle size	P/N	Price (yen)
XG-CN	I.D. 10.0x150mm		XGCN5P1150WN	130,000
	I.D. 10.0x250mm	h	XGCN5P1250WN	150,000
	I.D. 20.0x150mm		XGCN5P2150WN	240,000
	I.D. 20.0x250mm		XGCN5P2250WN	275,000

Prep guard column[Normal phase]

Packing name	Column Size	P/N	Price (yen)
XG-CN	(for I.D. 10.0mm)	XGCN580010WN	20,000
	I.D. 20.0x50mm	XGCN5P2050WN	105,000
The above goods are not cartridge type things.			
It becomes a packing guard column.			

At the time of the order, please appoint it for "a normal phase" or "a reverse phase" which it is.